





◆ Temperature and humidity control range diagram

◆ Inner tank

◆ 720FS

Name	Model	Temperature Range (°C)	Humidity Range (RH)	Capacity (L)	Interior Dimensions(mm) W×D×H	External Dimensions(mm) W×D×H	Power (kW)	Shelf (Standard)
Stability Chamber	Labonice-250FS	0 ~ 70	10 ~ 80%	250	600×500×830	780×880×1650	1.8	3/7
	Labonice-720FS	0 ~ 70	10 ~ 80%	700	975×575×1250	1280×880×1940	2.5	4/8
High-precision constant temperature test chamber	Labonice-250FT	0 ~ 70	N/A	250	600×500×830	780×880×1650	1.3	3/7
	Labonice-720FT	0 ~ 70	N/A	700	975×575×1250	1280×880×1940	2.0	4/8

All indicators on this color page were measured at an ambient temperature of 20~25°C.